

Small Animal Ultrasound Systems

LuciD™


D-CURVE™

Vi
PROBE™



See. Touch. Feel the difference.

4Vet Slim



LCD LED 10.4" Touchscreen



Quick and easy access to the key functions



Detailed, High-Resolution Scanning with Portability for Quick Emergency Screening, Orthopedic Surgeries and Companion Animal Examinations

Perfect images in all lighting conditions.

Veterinarian developed presets provide quick access to most common examination settings for multiple species and exam types.

Internal Battery for up to 4 h of constant work.

Biopsy needle echo enhancement function.

Starting at \$17,495 CAD.

Interchangeable Probes Available:

- Convex 2-5 MHz
- Convex 2-8 MHz
- Micro-Convex 4-9 MHz
- Linear 6-14 MHz 40 mm
- Linear Endorectal 4-9 MHz



BLUE



Vi-Probe™ - Virtual convex on linear probes,
 Virtual Phased Array on convex probes.

LuciD™ - Speckle reduction, contrast and edge
 enhancement system.

D-Curve™ - Grayscale curve perfectly matching
 the human eye perception.

Starting at \$19,495 CAD.

Interchangeable Probes Available:

- Convex R60 / 3.5 MHz
- Convex R50 / 5 MHz
- Linear L40 / 10 MHz
- Micro-convex R15 / 6.5 MHz
- P-probe R20 / 6.5 MHz
- Linear rectal probe L60 / 7 MHz
- OPU R11 / 6.5 MHz



LCD LED 12" Touchscreen



Ready to work in less than 30s



Internal battery allowing for 2.5h of constant work



SCANPad



Bright 15" Touchscreen Monitor



Intuitive and User-Friendly



Auto Image Optimization



Presets

OB, GYN, Cardiology,
Abdominal Settings

Modes

B, M, 2B, B/M, and 4B
scan modes

Starting at \$3665 CAD

CONSIDERING A NEW ULTRASOUND?

We'd love to be included for consideration because we believe (hand over heart) that we offer Canada's veterinarians the best in class, best value portable ultrasound scanners fully supported here in Canada.

Exceptional image resolution, features, support and models for all species and practices. Book your demo or trial now, talk to us!